

NMR Operator Training

TTKBL0004_EN

The procurement of the 700 MHz NMR spectrometer system was made in the framework of the GINOP-2.3.3-15-2016-00004 project supported by the EU and co-financed by the European Regional Development Fund.



Safety


- Personal safety
 - Cryogenics: do not lean on or push magnets
 - Cryoprobes: avoid contact with transfer line
 - Magnetic and RF hazards
- Instrument safety
 - Know the limits of instruments
 - Probe limits: avoid excessive long decoupling and long hard pulses or their equivalents
 - Be conservative
 - Double check pulse program and parameters for any non-standard new experiment
- Data Safety
 - Back up data promptly and regularly
 - Data processing or manipulation has no impact on the raw (FID) data
 - Do not change parameters after data are acquired

Two Insider Scoops

- Receiving efficiency
 - conceptually, it is similar to extinction coefficient in UV spectroscopy; it characterizes how efficient a unit magnetization can produce a signal by a given NMR receiver
 - NMR signal size is proportional to receiving efficiency
 - Receiving efficiency can be pre-calibrated as a function of 90° degree pulse length
 - receiving efficiency is the same for all nuclei of the same type (indifferent to chemical shifts) in the same sample
- Solvent signal offers a universal and robust concentration internal standard
 - Normalized NMR signal size is strictly proportional to concentration for a given sample, regardless how concentrated or dilute the sample is
 - Unit magnetization generates the same amount of total NMR response, which is indifferent to chemical shift or line-shape

Sample

- Rule #1: for Bruker NMR spectrometers, the NMR tube insert cannot exceed max depth (19mm or 20mm) from the center of the RF coil
 - Longer insert than recommended may present problems for the probe, as well as cause frictions during spinning
 - Varian is more flexible in allowing longer insert
- Rule #2: center of NMR sample should be as close as possible to the center of RF coil.
 - Normal sample needs to about 500 ul or slight more
 - Too much solvent is a waste!
 - Too little solvent may make shim difficult, but it does work!
- 10% deuterated solvent is sufficient for locking



Samples of smaller volumes

- Follow rule # 1 and then rule #2
- Shimming might be challenging due to air/glass and air/solution interfaces
- Consider Shigemi tubes
- Be careful with spinning
 - Non-spinning is recommended for volume ~ 300 ul or less




300ul

400ul

500ul


Sensitivity for smaller volumes

- Volume less than 300 ul may not offer additionally sensitivity improvement over that achieved by 300 ul, if the total amount of analyte is constant




Tune and match the probe

- Only higher fields (500, 600 and 800 MHz) in our facility need tuning
- Most of the time only proton channel requires tuning
- Drx500-2 with BBO needs special attention
 - Proton always needs tuning
 - BB (used for ^{13}C or ^{31}P etc) channel needs tuning, by first dialing the numbers to the pre-set values




Significance of tuning/matching


- Shorter 90 degree pulse
 - More efficient use of RF power
 - Protects transmitter
 - More uniform excitation in high power
- Better sensitivity
 - Reciprocity: if excitation is inefficient, then detection is equally inefficient
- Potentially quantitative:
 - The product of NMR signal size is inversely proportional to the 90 degree pulse length



Recognizing Bruker probe types



BBO probe on drx500-2



TXI probe

Lock

- Lock depends on shim: bad shim makes bad lock
 - Initialize shim by reading a set of good shims (i.e. rsh shims.txt)
 - Inheriting a shim set from previous users may present difficulties
 - Unusual samples (esp. small volumes) may need significant z1/z2 adjustments
- Use “lock_solvent” or “lock” command
 - The default (bruker) chemical shift may appear as dramatically changed if the spectrometer assumes another solvent
- Avoid excessive lock power
 - Lock signal may go up and down if lock power is too high due to saturation of deuterium signal
 - Apply sufficient lock power and gain so that lock does not drift to another resonance (this may happen by autolock if multiple deuterium signals exist)

Shim

- The goal of shim is to make the total magnetic field within the active volume homogeneous (preferably <1Hz).


Total magnetic field =

static field (superconductor) + cryoshim (factory set) + RT shim (user adjust)

- Shim can be done either manually or by gradient, which can be very efficient and consistent if done properly
- Sample spinning may improve shim
 - However, spinning-side band appear
- Recommendation:
 - Start from a known good shim set (by rsh on bruker or rts on varian).
 - Do not inherit shims from other users unless you know they're good
 - Non-spinning and higher order (spinning) shims should not change dramatically from sample to sample for most applications

Lock: lock gain


recommended



higher lock gain

Lower lock level
due to lower lock gain


not recommended



may easily lose lock;
change in lock level (during shimming) is less visible


Lock: avoid high lock power

Good lock



Lock power okay

Bad lock



Lock power too high


unstable lock and lower level

Evaluate shims

- Look for a sharp peak
 - No clear distortion
 - Full width at half height should be about 1 Hz or less for small molecules
 - Small (1% or smaller) or free of spinning side-bands
- Check if peak distortions are individual or universal
- Make sure that phasing is not causing peak distortions
- Maximize the lock level
 - Higher lock level => better shim
- Lock level does not drop significantly when spinning is turned off


Shim by line-shape

Plot made by G. Pearson, U. Iowa, 1991




Understanding NMR

- Modern NMR spectrum is an emission spectrum
- Equilibrium state
 - Magnetization is along +z axis
 - It is desired to have the largest +z magnetization prior to excitation
- Excitation by a RF pulse
 - A projection of magnetization is made on xy plane
 - It is desired to have the largest xy plane project for observation
- Observation
 - Precession of the projected xy-plane magnetization



RF pulses

- RF pulse manipulates spins
 - Important in excitation and decoupling
 - Defined by length, power and shape
- RF power is expressed in decibels
 - Bruker
 - Power range: typically 0db (high power) to 120db (low power)
 - Varian:
 - Coarse power: typically 60db (high) to 0db (low); 1 db increment; absolute
 - Fine power: 4095 (high) to 0 (low); default is 4095; relative
 - e.g. 54.5db can be roughly achieved through setting coarse power to 55 and fine power to 3854




RF pulse calibration

- Hard pulse (high power pulse) can be calibrated directly or indirectly
- For best calibrations, pulses need to be on resonance (know the chemical shift or resonance frequency!)
- Soft or shaped pulsed can be first calculated and then fine-tuned to optimum
 - Shapetool (by Bruker) or Pbox (by Varian) can be used for calculation and simulation
 - Be aware of possible minute phase shift (several degrees for soft pulses), which can be critical in water flip back or watergate

Proton pulse calibration


- Most hard (highest power) 90° pulses are typically from 5 us to 20 us.
- Direct observation for high power proton pulse calibration (or even for heteronuclei if sensitivity is sufficient)
 - 360° method (not quite sensitive to radiation damping or relaxation)
 - 180° method



NMR observables

- Chemical shifts
 - expressed in ppm
- Scalar couplings
 - expressed in Hz
 - 2D or nD bond correlations
- NOEs / relaxation / line-shapes
- Peak size
 - Potentially useful in quantitative analysis

From 1D to 2D



2D NMR

- Correlate resonances through bond or space
 - COSY: coupling
 - Magnitude mode recommended.
 - 1 mg or less will do
 - Minutes to a couple of hours
 - TOCSY: coupling network
 - ~ 70 ms mixing time
 - 1 mg or less will do
 - An hour or longer
 - NOESY / ROESY: distance / NOE
 - Mixing time ranging from less than 100 ms (proteins) to 500 ms (small molecules)
 - 1 mg or more
 - Hours or longer
 - HSQC/HMQC: proton correlation to X, typically through one-bond scalar couplings (two or three bond correlation possible)
 - 1mg or less will do
 - An hour or longer
 - HMBC: proton correlation to X, through multiple bond scalar couplings
 - 1 mg or more
 - Hours or longer

2D NMR

- Resolve overlapping peaks
 - Resolution is provided largely through the indirect dimension
 - No need to have highest resolution in the direct detected dimension
 - Limit direct acquisition time to 100ms or less if heteronuclear decoupling is turned on
 - Lower decoupling power if longer acquisition time is needed
 - Change in experimental conditions may help

2D NMR essentials: acquisition

- Proton tuning and matching
- Calibration of proton (90 degree) pulse length
 - Standard pulse lengths can be used if the solution is not highly ionic (< 50 mM NaCl equivalent)
 - All proton pulses are likely getting longer if the solution is ionic and/or the probe is not tuned
- Modest receiver gain
 - rg about half of what rga gives or less
- Carrier frequency (center of spectrum in Hz) and SW (sweep width) in both dimensions (avoid aliasing unless intended to)
- Number of scans (NS)
 - The pulse program recommends NS (a integer times 1, 2, 4, 8 or 16)
 - Needs some dummy scans, especially with decoupling / tocsy
- Number of increments in the indirect dimension (td1)
 - Larger td1 improves resolution in the indirect dimension
 - Rarely exceeds 512 (except occasionally in COSY)
- Detection method in the indirect dimension
 - Determined by the pulse program
 - Typically is either states (and/or TPPI) or echo-antiecho
- Acquisition time (aq) less than 100 ms with decoupling
- Modest gradients (cannot be more than the full power of 100% and typically less than 2 ms in duration)
- Go through the pulse program if you really care

2D processing

- Window functions
 - Allow FID approach zero at the end of the acquisition time
 - Sine bell functions with some shifts are recommended most of the time
- Zero filling
 - Typically double data points in each dimension
- Phasing
 - Indirect dimension zeroth and 1st order corrections are recommended in the pulprogram. If not, use 0 for both
 - Direct dimension first order phase is rarely more than 50 degrees. Zeroth order can be anywhere from 0 to 360 degrees
 - Phase in the 2D mode to best appearance
- Referencing
 - Can be done by picking a known resonance in the spectrum
 - Or referenced by (external) protons