



- D - can be heated directly
- ID - can be heated indirectly
- N - non heatable




test tube (D)




beaker (ID)




flat bottom flask (ID)




round bottom flask (D)




Erlenmeyer flask (ID)




distillation flask (D)




crystallization dish (ID)



volumetric flask (N)



measuring cylinder (N)




volumetric pipette (N)



graduated pipette (N)




burette (N)




weighing bottle (N)




watch glass (ID)




gas washing bottle (N)




U-tube (N)




funnel (ID)




analytical funnel (ID)




powder funnel (N)




Büchner funnel (N)




filter flask (N)




filter crucible (N)




filter crucible (N)




water jet pump (N)




separating funnel (N)



Liebig condenser (N)




bulb condenser (N)




spiral condenser (N)




desiccator(N)



evaporating dish (ID)




crucible with lid (D)




mortar with pestle (N)




double bosshead (clamp holder)




condenser clamp




extension clamp




flask clamp



burette clamp




pinchcock clamp Mohr pattern




filter ring




crucible tongs




pinchcock clamp Hoffmann pattern




wire gauze with asbestos (ceramic centre)




clay triangle




iron tripod




Bunsen support stand



support stand



Bunsen burner



washing bottle